






# FINOCHIETTO RETRACTORS



SKU	Size	Blade Width	Blade Height	Spread Width	Spread Height	Alloy
TL0071.1	Pediatric/Small	6mm	6mm	30mm	28mm	Stainless Steel
TL0072.1	Pediatric/Large	15mm	20mm	65mm	58mm	
TL0074.1	Adult/Small	35mm	40mm	145mm	135mm	Stainless Steel
TL0076.1	Adult/Medium	40mm	65mm	210mm	170mm	
TL0078.1	Adult/Large	60mm	65mm	210mm	190mm	