





POTTS SCISSORS

10MM SHORT FINE ANGLED BLADES



SKU	Shaft Angle	Handle	Total Length	Alloy
SL0715.1	25°	Round	6.75" (17cm)	Stainless Steel
SL1715.1			8" (20cm)	
SL0716.1	45°		6.75" (17cm)	
SL1716.1			8" (20cm)	
SL0717.1	60°		6.75" (17cm)	
SL1717.1			8" (20cm)	
SL0718.1	90°		6.75" (17cm)	
SL1718.1			8" (20cm)	
SL0719.1	125°		6.5" (16.5cm)	
SL1719.1			8" (20cm)	
SL0785.1	125° V-Neck	6.5" (16.5cm)		
SL0739.1		8" (20cm)		