



VATS SCISSORS



SKU	Jaw Type	Jaw Size	Shaft	Working Length	Total Length
SL8400.1S	Straight Razor Edge Blades w/Blunt Micro Tips	1cm	5mm	6.5" (16.5cm)	10.5" (26.5cm)
SL8402.1S	Curved Razor Edge Blades w/Blunt Tips	3cm	5mm	7.5" (19cm)	11.5" (29cm)
SL8404.1S	Curved Razor Edge Blades w/Blunt Tips	3.5cm	7.5mm	7.5" (19.5cm)	11.5" (29cm)