

VATS Foerster **Forceps**



- FL8436.12 Curved Left 12mm Oval Serrated jaws, 7mm Shaft, w/Ratchet, 9.5" (24cm) Working Length, Stainless Steel, 13.5" (34cm)
- FL8433.12 Curved Left 12mm Oval Serrated jaws, 7mm Shaft, w/Ratchet, 7.5" (19cm) Working Length, Stainless Steel, 11.5" (29cm)
- FL8437.12 Deeply Curved Left 12mm Oval Serrated jaws, 7mm Shaft, w/Ratchet, 9.5" (24cm) Working Length, Stainless Steel, 13" (33cm)

Thoracoscopic DeBakey **Clamp**



- AL8282.12 110mm Curved 1x2 DeBakey Atraumatic Tapered jaws, 10mm Shaft, w/Ratchet, 25.5cm Working Length, Stainless Steel, 14.5" (37cm)

VATS Node Grasping **Forceps**




- FL8429.12 Curved Left 11mm Oblong Serrated jaws, 5mm Shaft, w/Ratchet, 9" (23cm) Working Length, Stainless Steel, 13" (33cm)

VATS X-Heavy **Needle Holders**



- NL8410.15 Straight X-Heavy jaws w/TC Inserts, 10mm Shaft, w/Ratchet, 21cm Working Length, Stainless Steel, 13.5" (34cm)

MIS DeBakey **Forceps**



- FL8610.1 Straight 2x3 DeBakey jaws, 5mm Shaft, 10" (25cm) Working Length, Round Squeeze Handle, Stainless Steel, 14" (35.5cm)

VATS Mixer **Dissector**



- FL8418.1 Curved DeBakey 1x2 Atraumatic jaws, 5mm Shaft, 7.5" (19cm) Working Length, Stainless Steel, 11.5" (29cm)

VATS Right Angle **Dissector**



- FL8426.12 DeBakey/Cooley serrated jaws, 10mm shaft, 9.5" (24cm) Working Length, w/Ratchet, Stainless Steel, 13.5" (34cm)

Thoracoscopic **Scissors**




- SL8404.1S Curved 3.5cm Razor Edge blades w/Blunt tips, 7.5mm Shaft, 19.5cm Working Length, Stainless Steel, 11.5" (29cm)

Thoracoscopic **Knot Tier/Pusher**



- AL8450.12 5mm Shaft, 6.75" (17cm) Working Length, Round Squeeze Handle, Stainless Steel, 11.75" (30cm)

Wexler Metal **Sterilization Tray**



- DL2230 Single level w/Silicone Brackets, 10.5"x21"x2.5" (26.6x53.3x6.3cm)