



Designed in collaboration with Gustavo Bondulich, MD

VIPER FORCEPS

Features:

- Provide a secure grasp of lymph nodes during lymph nodes dissection
- Secure grasp of all kind of tissues like pleura, pericardium, diaphragm
- May be used as a dissector during VATS procedures
- The curved jaw improves visualization during VATS procedures
- Exclusive ratchet design allows tips to touch before the instrument locks
- May be used as a dissector during VATS procedures.



VATS

SKU	Jaw Type	Jaw Size	Shaft	Lock	Working Length	Total Length
FL8427.12	Curved Left	Tapered 5mm to 3mm	5mm	w/ Ratchet	9.5" (24cm)	14" (36cm)

