





JACOBSON MICRO SPIRAL NEEDLE HOLDERS



SKU	Suture Size	Tip Type	Lock	Handle	Length	Alloy
NL0710.21S	5-0 (and smaller)	Straight TC Coated Jaws	w/o Lock	Round	6.25" (16.5cm)	Titanium
NL0710.23S			w/Lock			
NL0712.21S			w/o Lock		7.25" (18.5cm)	
NL0712.23S			w/Lock			
NL0714.21S			w/o Lock		8.25" (21cm)	
NL0714.23S			w/Lock			
NL0716.21S			w/o Lock		9" (23cm)	
NL0716.23S						
NL0718.23S			w/Lock		12" (30.5cm)	Stainless Steel
NL0720.13S						
NL0720.21S			w/o Lock		10" (25cm)	Titanium
NL0720.23S			w/Lock			
NL0729.13S					7.25" (18.5cm)	Stainless Steel
NL0731.11S			w/o Lock			
NL0731.13S					8.25" (21cm)	
NL0731.13S-O			w/Lock			
NL0733.13S-O					9.25" (23cm)	