



# MIS BAYONET INSTRUMENTS



MIS Bayonet DeBakey Forceps

SKU	Tips	Shaft	Lock	Working Length	Total Length	Alloy
FL8706.1	Straight 2.8mm	5mm Bayonet Shaft	w/Lock	10.25" (26cm)	16.75" (42.5cm)	Stainless Steel



MIS Bayonet Needle Holder

SKU	Blades	Shaft	Lock	Working Length	Total Length	Alloy
NL8702.15	Curved Jaws w/TC Inserts	5mm Bayonet Shaft	w/Lock	10.25" (26cm)	16.75" (42.5cm)	Stainless Steel