




PEDIATRIC CLAMPS



Micro Coarctation Clamp

SKU	Jaw Type	Shanks Type	Total Length	Alloy
AL2179.1	Right Angled & Curved DeBakey Atraumatic jaws	Angled	4" (10.5cm)	Stainless Steel



Lambert-Kay Aorta Clamp

SKU	Jaw Type	Shanks Type	Total Length	Alloy
AL2126.1	Angled DeBakey Atraumatic jaws	Angled	8" (20cm)	Stainless Steel



Wexler Baby Vascular Clamp

SKU	Jaw Type	Shanks Type	Total Length	Alloy
AL2318.1-A	Profunda Cooley Atraumatic jaws	Deeply Angled	5.25" (14.5cm)	Stainless Steel