



# T-RIBBON RETRACTORS



SKU	End	Width	Alloy	Total Length
TL0350.1	1" (2.5cm) U Shaped	0.5" (1.1cm)	Malleable Stainless Steel	10.5" (26.5cm)
TL0351.1	2.26" (5.74cm) T Shaped	0.7" (1.78cm)		10.5" (26.7cm)
TL0352.1	1.2" (3cm) T Shaped	0.5" (1.1cm)		10.5" (26.5cm)
TL0353.1	2.75" (7cm) T Shaped	0.6" (1.52cm)		10.5" (26.7cm)
TL0354.1	1" (2.5cm) T Shaped	0.5" (1.25cm)		10.5" (26.5cm)
TL0355.1	3" (7.62cm) T Shaped	0.6" (1.52cm)		10.5" (26.7cm)
TL0356.1	1.5" (3.8cm) T Shaped	0.55" (1.4cm)		10.5" (26.5cm)
TL0357.1	3.27" (8.3cm) T Shaped	0.75" (1.9cm)		10.75" (27.3cm)