


VIPER FORCEPS

Features:

- Provide a secure grasp of lymph nodes during lymph nodes dissection
- May be used as a dissector during VATS procedures
- The curved jaw improves visualization during vats procedures
- Exclusive ratchet design allows tips to touch before the instrument locks
- Designed in collaboration with Gustavo Bondulich, MD



VATS

| SKU | Jaw Type | Jaw Size | Shaft | Lock | Working Length | Total Length |
|-----------|-------------|--------------------|-------|------------|----------------|--------------|
| FL8427.12 | Curved Left | Tapered 5mm to 3mm | 5mm | w/ Ratchet | 9.5" (24cm) | 14" [36cm] |

